Operating instructions

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Mode d'emploi

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
REF 224.40
REF 224.52

CE
0051

Pneumatic micromotor
Micromotore pneumatico
Micromoteur pneumatique
Micromotor neumático
Micromotor pneumático

TKD
TeKne Dental

ACCESSORIES AND SPARE PARTS / ACCESSORI E RICAMBI / ACCESSOIRES ET RÉCHANGES / ACCESORIOS Y RECAMBIOS / ACESSORIOS E PECAS SOBRESSELENTES

REF 001.89.01	O-rings / Guarnizioni / Joints toriques / Guarniciones / O-rings
REF 120.00	SO2000 Service Oil / Lubrificante / Lubrifiant / Lubricante / Óleo de serviço
REF 120.01	Nozzle / Beccuccio / Buse / Boquilla / Bocal
REF 120.02	Nozzle for rapid connection / Beccuccio per attacco rapido / Buse pour connexion rapide / Boquilla para conexión rápida / Bocal para conexão rápida.
REF 129.40	Air pressure gauge for 4-hole connection / Manometro per attacco 4 vie / Manomètre pour raccordement 4 voies / Manómetro para enlace 4 orificios / Manómetro de 4 vias
REF 316.40	GYROFLEX® rapid coupling for 4-hole connection / Raccordo rapido per connessione 4 vie / Raccord rapid pour raccordement 4 voies / Enlace rápido para conexión 4 orificios / Acoplamiento rápido para conexão 4 vias

DESCRIPTION

Pneumatic micromotor for professional dental use, compliant to ISO 13294 standard. The micromotor must be used with any straight or contra-angle handpiece with coupling meeting the requirements of ISO 3964 standard, with internal spray.

Classification: Type 2 according to ISO 13294.

Weight: 98 g.

Speed: variable from 5000 to 20000 rpm in both rotation modes.

Maximum torque: 3.5 Ncm.

Maximum noise level: 61 dBA.

This medical device meets the requirements of the European Directive 93/42 EEC (Class IIa).

CONTENTS

VORTICE® micromotor, a set of spare O-rings.

CONNECTION

Fixed connection (REF 224.40): 4-hole attachment connection meeting the requirements of ISO 9168 Type B standard.

Rapid connection (REF 224.52): for use with GYROFLEX® rapid couplings (REF 316.40). This version can also be used with the Multiflex® rapid couplings.

AIR SUPPLY

Recommended air pressure: 3.0 ± 0.1 bar (43.5 ± 1.5 psi).

Maximum air consumption: 63 NI/min.

Check air pressure by means of a gauge (REF 129.40) connected between the micromotor and the supply hose.

Supplied air must be dry and well-filtered, kept at constant pressure.

Impurity or residues in the compressed air may cause the instrument to malfunction.

SPEED

Speed may be increased or decreased by rotating the ring nut located on the body of the instrument (Fig. 1). In order to get a clockwise or anticlockwise rotation, turn the the ring nut to right or left respectively.

SPRAY

Maximum water pressure: 2.5 bar.

Water flow is higher than 50 ml/min at 2.0 bar.

The device is equipped with a non-return valve.

Water must be well-filtered.

LUBRICATION

In order for the instrument you have purchased to be helpful in your profession for a very long time, carefully repeat the following operation **before every sterilization cycle and, anyway, at least once a day.**

Lubricate by means of the special can (REF 120.00), after inserting the nozzle on it (REF 120.01), as shown on Fig. 3.

The SO2000 service oil must be inserted into instrument air input. Press once only.

In case of micromotor with rapid connection, insert and screw the special nozzle (REF 120.02) and lubricate as shown on Fig. 4.

After lubrication, run the instrument without handpiece for a few seconds .

REPLACEMENT OF O-RINGS

Should water or air leakage occur, replacement of O-rings on the coupling is recommended.

With the aid of a pin (Fig. 2), pull out the old O-rings and insert the new ones on the corresponding grooves (REF 001.89.01).

STERILIZATION

The micromotor can be sterilized in autoclave at 135°C / 2 bar (275°F/29 psi).

Before sterilization, remove residues and carefully clean the instrument.

The device is supplied "not sterile".

CLEANING AND DISINFECTION

The micromotor can be externally cleaned with a disinfectant containing: *phenol, gluteraldehyde up to 3%, ethyl alcohol.*

The chemical substance must not contain any chlorine, acid and solvent.

The micromotor cannot be immersed in any kind of liquid.

SERVICE

In case of malfunction and for any overhaul and repair work, please contact your usual supplier or directly to TKD.

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