**Operating instructions**  
**Istruzioni per l'uso**  
**Mode d'emploi**  
**Modo de empleo**  
**Instruções de uso**

REF 352.00

**TKD**

TeKne Dental

0051

Fig. 1



Fig. 2

### DESCRIPTION

Computer-controlled unit for use with TITANUS<sup>®</sup>E piezo-electric scaler.

The scaler can be used for:

- scaling prophylaxis procedures
- dry procedures (without spray cooling)
- endodontic procedures.

The unit:

- must be used only by trained and qualified personnel.
- is meant for continuous operation with intermittent loading
- does not produce electromagnetic interference
- is not suitable for use in the presence of a flammable anaesthetic mixture with air or with oxygen or nitrous oxide.

This medical device meets the requirements of the European Directive 93/42 EEC (Class IIa).

### CONTENTS

SELF-TITANUS<sup>®</sup> unit, power cord, pneumatic foot pedal, TITANUS<sup>®</sup>E piezo-electric scaler.

### SETTING UP

With reference to Fig. 1, connect the following parts:

- a) Water input. Connect with 4 mm I.D. tubing. Use well-filtered water or disinfectant. In any case, connecting the supplied water filter is recommended.
- b) Foot control input.
- c) Mains input.

Additionally, the following parts are shown:

- d) Protection fuses (2 x F 500 mA).
- e) Power switch

### USE OF THE HANDPIECE

For connection, use and maintenance of the TITANUS<sup>®</sup>E scaler, refer to corresponding operating instructions.

To avoid any risk of injury, the insert must always be removed before resting the handpiece on the unit's holder.

### POWER ADJUSTMENT

- ☐ Power can be adjusted by means of the two upper push-buttons (Fig. 2). Setting ranges from 1 (minimum power) to 10 (maximum power). Values of settings are shown on the LCD as long as the handpiece is connected.

### WATER ADJUSTMENT

- ☐ Water flow can be adjusted by means of the two lower push-buttons (Fig. 2). Setting ranges from 0 (flow closed) to 10 (maximum flow). Value of setting is shown on the display. Water flow is affected by water input pressure. Cooling spray flow rate is higher than 50 ml/min. (@ 0.5 bar and maximum flow setting). Values of settings are shown on the LCD as long as the handpiece is connected.

### ACTIVATION

Vibrations of scaler are activated by the pneumatic on/off foot control.

### ENDO MODE

Decreasing power down to 1 and keeping the ☐ push-button pressed for about two seconds will let you enter the Endo mode: now power is reduced to about 50% of normal Scaler level. In order to re-enter Scaler mode, increase power up to 10 and keep the ☐ push-button pressed for about two seconds.

### CLEANING AND DISINFECTION

Cleaning of control box: by means of ordinary detergent products.  
Cleaning of the silicone hose: by means of a disinfectant based on alcohols like *ethanol*, *methanol* or a product containing up to 3% of *glutaraldehydes*.

### GUARANTEE

Either the unit or the handpiece carry a 12-month guarantee against all defects of material or faulty construction.

### SERVICE

In case of malfunction or failure of the integrity of the enclosure, and for any overhaul and repair work, please contact your usual supplier or directly to TKD.

### REFUSE DISPOSAL

- ♻️ This device needs to be recycled. Electrical and electronic equipment may contain dangerous substances which constitute health and environmental hazards. The user must return the equipment to its dealer or establish direct contact with an approved body able to process and derive value from this type of equipment (European Directive 2002/96/EC).

### TECHNICAL DATA

Power supply: 230 V~ (alternating current), 50/60 Hz, 21 VA

Operating frequency: 28000 ÷ 31000 Hz

Conformity: CEI EN 60601-1, CEI EN 60601-1-2

- ♻️ Degree of protection against direct and indirect contacts: Class I, type B applied part according to CEI EN 60601-1.

Degree of protection against ingress of water: ordinary IPX0 equipment.

Maximum water input pressure: 4 bar (58 psi).

Sizes: 150 x 134 x 65 mm

Weight without handpiece: 1100 g

The scaler is meant for intermittent use with the following working times:

- with active spray cooling and at minimum power: maximum 30 min. ON, minimum 10 min. OFF
- without spray cooling and at minimum power: maximum 10 min. ON, minimum 10 min. OFF
- with active spray cooling and at maximum power: maximum 8 min. ON, minimum 10 min. OFF
- without spray cooling and at maximum power: maximum 6 min. ON, minimum 10 min. OFF

Working temperature: 0 ÷ 40 °C (32 ÷ 104 °F)  
Storage temperature: -20 ÷ 100 °C (-4 ÷ 212 °F)  
Working humidity: 30 ÷ 90 %  
Storage humidity: 10 ÷ 90 %

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