



# GYROFLEX®LED

REF 316.62  
REF 316.62.09

Rapid couplings  
Raccordi rapidi  
Schnellkupplungen  
Raccords rapides  
Enlaces rápidos  
Acoplamentos rápidos  
Ταχυσύνδεσμος

Operating instructions  
Istruzioni per l'uso  
Bedienungsanleitung  
Mode d'emploi  
Modo de empleo  
Instruções de uso  
Οδηγίες χρήσεως

ACCESSORIES AND SPARE PARTS / ACCESSORI E RICAMBI / ZUBEHÖR UND ERSATZTEILE / ACCESSOIRES ET RÉCHANGES / ACCESORIOS Y RECAMBIOS / ACCESÓRIOS E PEÇAS SOBRESSELENTES / ΕΞΑΡΤΗΜΑΤΑ ΚΑΙ ΑΝΤΑΛΛΑΚΤΙΚΑ

REF 026.00.00  
REF 001.29.01

LED  
Green O-rings / Guarnizioni verdi / Grüne O-Ringe / Joints toriques vertes / Guarniciones verdes / Orings verdes / Πράσινοι ελαστικοί δακτύλιοι (10 pcs.)

REF 001.29.02

Black O-rings / Guarnizioni nere / Schwarze O-Ringe / Joints toriques noirs / Guarniciones negras / Orings pretos / Μαύροι ελαστικοί δακτύλιοι (10 pcs.)

## ENGLISH

### DESCRIPTION

Rapid couplings for use with all quick-release TKD turbines or with any other dental instrument with Multiflex® connection (Fig. 1).

Couplings are equipped with a non-return valve for water.

Maximum operating pressure of air and water: 6 bar (87 psi).

The rapid coupling is provided with an integrated LED (3.5 ± 0.2 Vdc; 60 mA maximum) for illumination.

This medical device meets the requirements of the European Directives 93/42/EEC (Class I) and 2007/47/EC.

### CONTENTS

Rapid coupling, tool for coupling installation, set of spare O-rings.

### CONNECTION TO THE SUPPLY HOSE

Screw the coupling onto the flexible supply hose. The connector meets ISO 9168 Standard.

For a perfect locking, use and insert the supplied tool on the corresponding grooves (Fig. 2).

After that, the coupling does not need to be removed from the hose (Fig. 3).

Couplings must be connected to the power supply unit meeting the requirements of CEI EN 60601-1 Standard.

Additionally, the power supply unit must have a minimum insulation strength of 500 V between the electric circuits and earth and a minimum insulation strength of 4000 V between the electric circuits and mains.

### SPRAY REGULATION

Rotate the spray ring on the coupling in order to regulate the fraction of water (Fig. 5). Different stop positions allow a precise regulation of the water volume. Turn the ring clockwise to reduce the water volume. Turn the ring counterclockwise to increase the water volume.

### LED POSITIONING

The LED must be supplied with the right polarity. If the LED does not turn on: unscrew the protection cap (a), pull out the LED (b), rotate it by 180° so as to reverse polarity on its contacts and insert it again on its housing (c) (Fig. 4). Screw cap again.

### O.R. MAINTENANCE

In order to increase coupling seal and smoothness, we recommend that you regularly rub the O-rings with a small quantity of non-toxic grease. The grease can also be taken from the turbine grease lubricator (LubriONCE®)

### O.R. REPLACEMENT

Should water or air leakage occur, replacement of O-rings on the coupling is recommended.

With the aid of a pin, pull out the old O-rings and insert the new ones on the corresponding grooves (Fig. 6).

### CLEANING AND DISINFECTION

Couplings can be externally disinfected with a non-chlorinated germicidal cleanser. Additionally, the cleanser must not contain any glutaraldehyde, acid or solvent.

Rapid couplings cannot be sterilized.

Multiflex® is a registered trademark of Kaltenbach & Voigt GmbH, Germany.

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