### DESCRIPTION

High-torque ceramic ball-bearing turbine for professional dental use, compliant to EN ISO 14457 standard.  
 Weight without coupling: 38 g.  
 Maximum speed of rotation: 345000 rpm.  
 Output illumination: 25000 lux.  
 Maximum noise level: 58 dBA (@3 bar).

This medical device meets the requirements of the European Directives 93/42/EEC (Class IIa) and 2007/47/EC.

### CONTENTS

BRAVIA®XL turbine, a LubriONCE® grease lubricator, a wire cleaner.

### CONNECTION

The BRAVIA®XL turbine must be used together with the GYROFLEX® (REF 316.60) or GYROFLEX®LED (REF 316.62 or REF 316.62.09) rapid couplings (Fig. 1). These turbines can also be used together with the MULTIflex® rapid couplings.

### AIR PRESSURE

**Recommended air pressure: 2.8 ± 0.2 bar (40.5 ± 3.0 psi).**  
 Air consumption: 43 Nl/min (@3 bar).  
 Check air pressure by means of a gauge (REF 129.40) connected between the turbine and the supply hose (Fig. 2).  
 Supplied air must be perfectly dry and well-filtered, kept at constant pressure.

### SPRAY

Three-jet type.  
 Spray air pressure: 1.0 ÷ 4.0 bar (14.5 ÷ 58.0 psi).  
 Water pressure: 0.8 ÷ 2.5 bar (11.6 ÷ 36.3 psi). Water must be well-filtered.  
**The use of unfiltered hard water will lead to early blockage of the tubings and spray diffusers.**  
 To clean the spray ducts, use the wire cleaner supplied, as shown on Fig. 4.

### LIGHTING

The BRAVIA®XL turbine is equipped with a highly-efficient light conductor. The light source is located in the GYROFLEX® or GYROFLEX®LED coupling.

### BLUR TYPE

Shank diameter: Ø1.60 mm, Type 3 according to ISO 1797-1.  
 Maximum length: 25 mm, code 4 or 5 or 6 (6<sup>th</sup> digit) according to ISO 6360-1.  
 Minimum fitting length of shank: 11 mm.

### BLUR CHANGING

The turbines are equipped with a sophisticated push-button chuck: press the button with the thumb and subsequently insert or extract the bur (Fig. 3).

### LUBRICATION BY GREASE

The turbine must be lubricated with the supplied LubriONCE® grease lubricator daily. At the very beginning, remove the cap and turn the rear housing clockwise until you can see grease coming out from the rod.

### LUBRICATION BY OIL

As an option to grease, it is possible to lubricate instrument by means of the SO2000 service oil (REF 120.00) carrying out the following operations:

- 1) Disconnect turbine from the rapid coupling;
- 2) Insert and screw the special nozzle (REF 120.02) onto the lubricant can;
- 3) Holding turbine by means of a cloth, fully insert the nozzle into the turbine (Fig. 6). Lubricate by pushing once only for one second.

In order for the instrument you have purchased to be helpful in your profession for a very long time, carefully repeat the described operations **before every sterilization cycle and, anyway, at least TWICE A DAY.**

### PRELIMINARY PREPARATION

Maintenance processes are described according to ISO 17664 standard. Before starting a process:

- Observe suitable personal protective measures against risk of infection.
- Remove all residual blood, cement or composite.
- Remove the bur from the instrument.
- Clean the medical device immediately after treatment of each patient. Do not immerse the medical device into disinfectants solutions and do not use ultrasonic cleaning units.

### CLEANING

Manual external: Carefully clean the surfaces using a clean cloth soaked in a suitable substance. Use also a medium toothbrush and brush off under flowing hot tap water.

### DISINFECTION

Manual external: Carefully clean the surfaces using a clean cloth soave in a suitable substance (i.e. O-phenilphenol or alcohols based on ethanol). Products not recommended: disinfectants containing benzalkonium chloride, acetone or glutaraldehyde.

### DRYING

Manual: In general use dry and clean air complying to ISO 7494-2 standard.  
 Automated: The drying phase is normally part of the cleaning program of the thermodesinfector.

### STERILIZATION

Sterilize with a steam autoclave with fractioned initial vacuum phase, class B cycle according to the EN 13060 standard. Sterilize at the nominal temperature of 134 °C (273 °F) for at least 3 min. Use only demineralised or distilled water. The device is supplied "not sterile". The LubriONCE grease lubricator and GYROFLEX® rapid couplings cannot be sterilized.

### STORAGE

The processed device must be kept protected from dust and microbial contamination. Storage temperature: -10 ÷ 70 °C. Storage humidity: 10 ÷ 90 %.

### GUARANTEE

The device carry a 24-month guarantee against all defects of construction.

### SERVICE

In case of irregular noise, head vibrations, bur swinging or any other malfunction of the instrument, and for any overhaul and repair work, please contact your usual supplier, an authorized service centre, or directly to TKD.

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