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ENGLISH

DESCRIPTION

The sophisticated USC60 electronic board (REF 387.00) allows to control operating power of any piezo-electric scaler working at any frequency comprised in the 27 ÷ 32 KHz band.

The board automatically adjusts operating frequency of vibrations to the real resonance load, which changes with the quantity of water delivered and the type of insert used.

Once wired to the USC60 electronic board, either the insert and the scaler handle are galvanically insulated, thus assuring absolute safety to the user.

CONTENTS

USC60 electronic board, two plastic screws for board assembly.

ASSEMBLY

Install and fix USC60 heat sink on the frame of the Unit using the two supplied screws and leaving the board sufficiently far away (at least 10 mm) from other components and/or the Unit frame.

Connection with TITANUS®E scaler (REF 400.00.0x or REF 410.00.0x) or TITANUS®S (REF 450.00.0x or REF 460.00.0x): connect red and black wires coming from the supply silicone hose (REF 330.1x, REF 332.1x or REF 333.1x) to the USC60 connector, as shown on attached wiring diagrams.

Board is activated by closing the START-STOP switch. Board can also be activated by an external voltage V_{st} (TTL signal, active low).

The board cannot be activated if the handpiece is not connected to the supply silicone hose.

A 24 V solenoid valve for water flow control can be connected as shown on the attached wiring diagrams.

 To protect the board it is necessary to connect an external fuse (F 2 A) as shown on the attached wiring diagrams.

POWER ADJUSTMENT

Power can be adjusted by means of a electric potentiometer (4,7 K Ω linear) or by an external continuous voltage V_{pw} (0 ÷ 5 Vdc). Please see attached wiring diagrams.

If the scaler is to be used for Endodontic procedures, output power can be limited to 50 %: this is accomplished by closing a switch from the connector no. 4 to the common ground (connector no. 6).

SERVICE

In case of malfunction and for any overhaul and repair work, please contact an authorized technician or apply directly to TKD.

REFUSE DISPOSAL

 This device needs to be recycled. Electrical and electronic equipment may contain dangerous substances which constitute health and environmental hazards. The user must return the equipment to its dealer or establish direct contact with an approved body able to process and derive value from this type of equipment (European Directive 2002/96/EEC).

DECLARATION OF RESPONSIBILITY

The manufacturer will not be held responsible in case of:

- use on a electrical system not meeting the requirements of CEI EN 60601-1 Standard
- operations or repairs carried out by technicians not authorized by the manufacturer

USC60 SPECIFICATIONS

The device:

- does not produce electromagnetic interference and meets the requirements of CEI EN 60601-1-2 Standard
- is not suitable for use in the presence of a flammable anaesthetic mixture with air or with oxygen or nitrous oxide
- must be connected to a power transformer having a minimum insulation strength of 500 V between secondary winding and earth and a minimum insulation strength of 1500 V between secondary winding and mains

Standard Conformity:	CEI EN 60601-1, CEI EN 60601-1-2, IEC 80601-2-60
Protection against electrical shock:	type B applied part according to IEC 80601-2-60 
Power supply:	24 Vac \pm 10 % (50/60 Hz) or 32 Vdc \pm 15 %
Maximum consumption:	24 VA at 24 Vac, 1 A at 32 Vdc
External protection fuse:	F 2 A
Operating frequency:	27 ÷ 32 KHz
Maximum Output voltage:	125 V_{RMS}
Maximum Output power:	9 W
Absolute maximum Input voltage:	40 Vdc
Operating Mode:	Continuos
Sizes:	60 x 54 x 30 mm
Working temperature:	0 ÷ 40 °C (32 ÷ 104 °F)
Storage temperature:	-10 ÷ 70 °C (14 ÷ 158 °F)
Working humidity:	30 ÷ 90 %
Storage humidity:	10 ÷ 90 %
Guarantee:	2 years

Protection: main oscillator is shut off in the following events:

- 1) Breaking of one supply wire
- 2) Short circuit between handpiece supply wires
- 3) General malfunction of the electronic control

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