TeKne Dental

APPLICATION NOTES

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January 2005

Gennaio 2005

LC2 ELECTRONIC BOARD

SCHEDA ELETTRONICA LC2

TeKne Dental s.r.l.

Via del Pescinale 77 - 50041 Calenzano (FI) - ITALY
info@teknedental.com www.teknedental.com

DESCRIPTION

The LC2 electronic board (REF 381.00) provides a universal power supply for two bulbs of all types of fibre-optic instruments. Each bulb can be turned on by an external electric switch. In order to drive the board by the instrument's air supply, the switch can also be pneumatic. Activation of LAMP 1 takes priority over the other.

Input voltage has a wide range and can be either alternative or direct, whereas output voltage (lamp power supply) is direct and can be adjusted by means of a trimmer.

By removing one jumper bulb turn-off time can be delayed up to 10 seconds from the corresponding switch aperture (Fig. 2).

ASSEMBLY

The board can be fixed to the desired location by means of the four supplied nylon supports. It is only necessary to make four holes as shown on Fig. 1.

CONNECTION

To connect board properly, please follows wire connections shown on Fig. 2 or 3:

- On Fig. 2 the wiring diagram for both bulbs are shown;
- On Fig. 3 the wiring diagrams for power supply of one bulb only are shown, the former with the activation switch on the supply circuit and the latter on the output circuit respectively; in both cases delayed turn-off is disabled.

To protect the board, it is necessary to connect an external fuse (F 250 mA) as shown.

Board must be connected to a power transformer having a minimum insulation strenght of 500 V between secondary winding and earth and a minimum insulation strenght of 4000 V between secondary winding and mains.

The manufacturer will not be held responsible in case of use on a electric connection not meeting the requirements of IEC 601-1 Standard.

SPECIFICATIONS

- Input voltage: 14 ÷ 26 Vac or 19 ÷ 36 Vdc.
- Output voltage: 2.6 ÷ 3.8 Vdc.
- Maximum output current: 1 A
- Maximum consumption: 3.1 W (Bulb) + 1.1 W (Board).
- Lamp turn-off delay time: 5 s or 10 s selectable by a jumper.
- External protection fuse: F 250 mA
- Circuit activation by means of low-current switches.
- Board sizes: 72 x 46 mm
- Working conditions: Temperature: 0 ÷ 40 °C (32 ÷ 104 °F); Humidity: 30 ÷ 90 %
- Storage conditions: Temperature: -20 ÷ 100 °C (-4 ÷ 212 °F); Humidity: 10 ÷ 90 %
- Guarantee: 2 years.

STANDARDS REFERENCE

The board does not produce electromagnetic interference and meets the requirements of IEC 601-1 and IEC 601-1-2 Standards.

SERVICE

In case of malfunction and for any overhaul and repair work, please contact an authorized technician or directly to TKD.

The manufacturer will not be held responsible in case of operations or repairs carried out by technicians not authorized by the manufacturer.

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