
Contra-angle Handpiece
Manipolo Contrangolo
Winkelstück
Pièce à main contre-angle
Pieza de mano contra-ángulo

REF 650.00
 REF 652.00

Operating instructions
Istruzioni per l'uso
Gebrauchsanleitung
Mode d'emploi
Modo de empleo

DESCRIPTION

Contra-angle handpiece with internal spray, direct 1:1 ratio (blue) and push-button grip. Handpiece compliant to EN ISO 14457 standard.
 Expected use: device intended for professional use in Dentistry for prophylaxis, restoration and endodontic works.
 This medical device meets the requirements of the European Directives 93/42/EEC (Class IIa) and 2007/47/EC.
 Maximum rotation speed: 40000 rpm.
 Weight: 63 g.

CONTENTS

QUARK®CA or QUARK®CAL handpiece, a wire cleaner.

CONNECTION

The handpiece can be used with any pneumatic or electric micromotor with internal spray and coupling compliant to ISO 3964 standard.
 Insert the handpiece onto micromotor's coupling and push it fully. In case of connection to a micromotor with light, rotate the handpiece until it clicks into place: in this way the fiber-optics rod will be aligned to the micromotor's light source.

Do not connect the handpiece to the micromotor while it is rotating.

SPRAY

One-jet type.
 Spray air pressure: 1.0 ÷ 3.0 bar (14.5 ÷ 43.5 psi). Air consumption is higher than 1.5 Nl/min (@2 bar).
 Spray water pressure: 0.8 ÷ 2.0 bar (11.6 ÷ 29.0 psi). Water consumption is higher than 50 ml/min (@2 bar). Water must be well-filtered.
 The device is equipped with a non-return valve for water.

BUR TYPE

Shank diameter: Ø2.35 mm, type 1 according to ISO 1797 standard.
 Maximum length: 22 mm, code 4 according to ISO 6360-1 standard.
 Minimum fitting length of shank: 12 mm.
 Bur is not part of the device. It is recommended to employ biocompatible burs compliant to ISO 10993-1 Standard.

REMOVAL AND LOCKING OF THE BUR

To lock a bur: push the button and insert the bur.
 To remove the bur: push the button and pull out the bur (Fig.1).

When fitting on a bur make sure it is fully inserted and correctly in place.

The device must not be started without a bur inserted into the chuck.

To avoid overheating of the push-button, which could cause burns, the button should not be pressed while the turbine is running. Internal tissues (tongue, cheeks, etc..) must then be protected by a retractor or a dental mirror.

PRELIMINARY PREPARATION

After each treatment, within a maximum of 30 minutes, clean, disinfect and lubricate the instrument, then sterilize it. Observing this procedure eliminates any blood, saliva or saline solution residues and prevents the transmission system from being blocked.

Maintenance processes are described according to ISO 17664 standard.

Before starting a process:

- Observe suitable personal protective measures against risk of infection.
- Remove all residual blood, cement or composite.
- Remove the bur from the instrument.
- Clean the medical device immediately after treatment of each patient.

Do not immerse the medical device into disinfectants solutions and do not use ultrasonic cleaning units.
 Do not immerse in physiological liquid (NaCl).

CLEANING

Manual external: Carefully clean the surfaces using a clean cloth soaked in a suitable substance. Use also a medium-hardness toothbrush to brush off residues.
 To remove any incrustations in the spray nozzle, use the supplied wire cleaner, as shown on Fig. 2.

Manual internal: Using a suitable spray product, follow corresponding instructions and spray the substance to the internal of the medical device. Immediately after internal cleaning, carry out the disinfection process and/or the sterilization process.

Automated: Use a thermidisinfector complying with ISO 15883 standard which is operated with alkaline cleaning agents at a pH of maximum 10.

DISINFECTION

Manual external: Carefully clean the surfaces using a clean cloth soave in a suitable substance (i.e. O-phenilphenol or alcohols based on ethanol). Products not recommended: disinfectants containing benzalkonium chloride, acetone or glutaraldehyde.

Manual internal: Using a suitable spray product, follow corresponding instructions and spray the substance to the internal of the medical device. Immediately after internal disinfection, carry out the OIL lubrication process and/or the sterilization process.

Automated: Use a thermidisinfector complying with ISO 15883 standard which is operated with alkaline disinfectant agents at a pH of maximum 10.

DRYING

Manual: In general use dry and clean air complying to ISO 7494-2 standard.

Automated: The drying phase is normally part of the cleaning program of the thermidisinfector.

LUBRICATION

It is possible to lubricate instrument by means of the SO2000 service oil (REF 120.00) carrying out the following operations:

- 1) Insert and screw the special nozzle (REF 120.03) onto the lubricant can;
- 2) Detach handpiece from the micromotor;
- 3) Holding handpiece by means of a cloth, fully insert the nozzle into the handpiece (Fig. 3). Lubricate by pushing once only for one second.

In order for the instrument you have purchased to be helpful in your profession for a very long time, carefully repeat the described operations **before every sterilization cycle and, anyway, at least TWICE A DAY.**

STERILIZATION

The device is supplied "not sterile". Sterilize it before first use.

Sterilize with a steam autoclave with fractioned initial vacuum phase, class B cycle according to the EN 13060 standard. Sterilize at the nominal temperature of 134 °C (273 °F) for at least 3 minutes. Use only demineralised or distilled water.

The device can be sterilized for at least 1000 times.

The quality of the sterilization depends very much on the cleanliness of the device. Only perfectly clean devices should be sterilized.

After cleaning, disinfecting, lubricating and sterilizing the device, before using it on the patient, insert a bur and run it at low speed for about 10 seconds, in order for the excessive lubricant to spread and exhaust.

STORAGE

The processed device must be kept protected from dust and microbial contamination.

Storage temperature: -10 ÷ 70 °C (14 ÷ 158 °F)

Storage humidity: 10 ÷ 85 %

Athmosferic pressure: 50 ÷ 106 KPa

GUARANTEE

The device is guaranteed against all defects of manufacture for a period of 24 months. Guarantee validity is subject to specific conditions.

SERVICE

In case of irregular noise, head vibrations, bur swinging or any other malfunction of the instrument, and for any overhaul and repair work, please contact your usual supplier, an authorized service centre, or apply directly to TKD.

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