

BUILT-IN SET Electric Micromotor Micromotore elettrico

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DESCRIPTION

Brushless 3-phase electric micromotors for professional dental use, compliant to IEC 80601-2-60 Standard.
Intended use: micromotors must be employed with any straight or contra-angle handpiece with internal spray and connection compliant to ISO 3964 Standard.
The micromotor can be used for prophylaxis, general dentistry and endodontics works. The NUCLEUS®LED model is equipped with a LED light source and can be used with any straight or contra-angle handpiece provided with fibre-optics.
Furthermore, the used handpiece must be biocompatible according to ISO 10993-1 Standard.

This medical device meets the requirements of the European Directives 93/42/EEC (Class IIa) and 2007/47/EC.

The device:
- must be used by competent and qualified personnel only;
- does not produce electromagnetic interference;
- is not suitable for use in the presence of a flammable anaesthetic mixture with air or with oxygen or nitrous oxide;
- is not liquid-proof (IP 20).

CONTENTS

NUCLEUS®LED micromotor with a set of spare O-rings, silicone hose with rotating connector, BMC60 electronic board.

CONNECTION

The NUCLEUS®LED (REF 612.00) micromotor must be used together with the specific silicone hose with rotating connector (REF 324.9x.01) and with the specific BMC60 (REF 397.00) electronic board.
Insert hose connector (Fig. 1) into the micromotor rear part and screw the metal nut. Once connected, the micromotor can be rotated through $\pm 270^\circ$ (Fig. 2).

Fit the straight or contra-angle handpiece onto the micromotor's coupling and insert it fully.

Never insert the handpiece when the micromotor is rotating.

INSTALLATION

To correctly install and connect the micromotor and the silicone hose to the specific power supply unit, please refer to the corresponding Application Notes of the BMC60 electronic board.
The power supply unit must meet the requirements of the CEI EN 60601-1 standard.

On final installation, check dielectric strength and leakage currents according to CEI EN 60601-1 standard.

On installation, micromotor must be sterilized before its first use.

COOLING

Micromotor must be cooled by means of compressed air coming from the dental unit. The minimum recommended air flow is 20 Nl/min while the minimum recommended air pressure is 2.5 bar (36 psi).
Maximum cooling air pressure: 5.0 bar (72 psi)

SPRAY

Maximum spray air pressure: 4.0 bar (58 psi). At 2.0 bar air flow rate is higher than 1.5 l/min.
Maximum spray water pressure: 4.0 bar (58 psi). At 2.5 bar water flow rate is higher than 50 ml/min.
The device is equipped with a non-return O-ring for spray water.
Water must be wellfiltered.

CLEANING AND DISINFECTION

After each treatment, clean and disinfect the device immediately. Observing this procedure can easily eliminate any residues of blood, saliva, cooling spray or other.
The exterior body of the micromotor can be cleaned with a disinfectant containing:
- glutaraldehyde up to 2%
- O-phenylphenol
- alcohols based on ethanol
Disinfectants containing chlorine, acids or solvents are not recommended.
The device cannot be immersed in any kind of liquid.

Micromotor does not need any lubrication.

STERILIZATION

Sterilize with a steam autoclave with fractioned initial vacuum phase, class B cycle according to the EN 13060 standard.
Sterilize at the nominal temperature of 134 °C (273 °F) for at least 3 min.
Use only demineralised or distilled water.
The device is supplied "not sterile".

O.R. REPLACEMENT

Should water or air leakage occur, replacement of O-rings on the coupling is recommended: with the aid of a pin (Fig.3), pull out the old O-rings and insert the new ones (REF 001.89.01 and REF 001.98.01) on the corresponding grooves.

SERVICE

In case of malfunction of the device, do not try to disassemble micromotor. For any check-up, overhaul or repair work, please contact your supplier, an authorized centre or directly TKD.

GUARANTEE

The device is guaranteed against all defects of manufacture for a period of 24 months. Guarantee validity is subject to specific conditions.

REFUSE DISPOSAL

This device needs to be recycled. Electrical and electronic equipment may contain dangerous substances which constitute health and environmental hazards. The user must return the equipment to its dealer or establish direct contact with an approved body able to process and derive value from this type of equipment (European Directive 2002/96/EC).

TECHNICAL SPECIFICATIONS

Standard conformity: EN ISO 14457: 1997, IEC 80601-2-60: 2012, CEI EN 60601-1: 2007, CEI EN 60601-1-2, European Directives 93/42/EEC and 2007/47/EC.

Classification: class 2, according to CEI EN 60601-1 standard

Protection against electrical shock: type B applied part according to IEC 80601-2-60 standard.

Maximum rotation speed: 40000 rpm
Minimum rotation speed: 100 rpm
Direction of rotation: clockwise or anti-clockwise
Maximum torque: 3.5 Ncm
Maximum peak current: 7 A (5 s)
Maximum noise level: 51 dBA (@45 cm)
Cooling air consumption: 20 Nl/min @ 2.5 bar
Sizes: $\varnothing 21.8 \times 73.6$ mm
Weight: 81 g
Working temperature: 0 ÷ 40 °C (32 ÷ 104 °F)
Working humidity: 30 ÷ 90 %
Storage temperature: -10 ÷ 70 °C (14 ÷ 158 °F)
Storage humidity: 10 ÷ 90 %
Light source: LED

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